
**Instructions for Use / Mode d'emploi /
Instrucciones de uso**

Background: The DynaCleft Nasal Elevation System (DNES) provides gentle support to the soft tissue of the nose. This support results in lengthening the columella and improving the preoperative shape of the nose. The DNES is typically used for 10-12 weeks or as indicated by the referring medical professional. All components and packaging do not contain natural latex rubber.

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**Safety Precautions / Consignes de sécurité /
Precauciones de seguridad**

1. Use a fresh DynaCleft Nasal Elevator System daily.
2. Do not use if the packaging has been opened or tampered with.
3. Never leave the device loosely attached to the baby.
4. The suspension tape should always be tight enough for the hook to stay in the nostril – replace the adhesive tether if the hook starts to fall off or dislodge.
5. Always use the device with a skin barrier tape affixed to the forehead.
6. Ensure baby sleeps on his/her back.
7. Remove the device when baby is in a prone (face-down) position.
8. Ask your doctor immediately about applying arm restraints when baby starts to grab, or roll from their back onto their stomach.
9. Arm restraints and/or mittens may be suggested to prevent baby from interfering with device location and tension. Consult your doctor.
10. Never leave other children alone or unsupervised with the baby.

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EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Made in Canada / Fabriqué au Canada / Hecho en Canadá

The DynaCleft Nasal Elevation System is used under medical supervision and is indicated for use in:

- Pre-surgical treatment of the nose prior to the first surgical repair of the cleft lip. This non-invasive device typically is used for as long as possible before the infant's first cleft repair surgery. The DNES has been designed to be used with the DynaCleft Uni-Lateral and Bi-lateral closure systems. Consult your medical professional if this is applicable.
- It may also be recommended by the medical professional as post-surgical stabilization and support after the first lip repair surgery.

Do Not Reuse: **The DNES system is intended to be replaced daily.**

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1. Repeated application and removal of any adhesive product can irritate sensitive skin.
2. Using more than the recommended tension may result in skin damage or adhesive failure.
3. Single use product not to be reused due to the effect on performance, cross contamination and risk of infection. Reuse can degrade adhesion strength and reduce effectiveness.

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What you will need

- DynaCleft Nasal Elevation System: hook, adhesive tether, paper template
- Scissors
- Recommended: A skin barrier dressing or tape (purchased separately). Consult your medical professional for advice on appropriate products.

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
A) Use the paper template provided to custom size the adhesive tether to the baby's nose-to-forehead dimensions. Thread tail of template through opening in hook, as shown.

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
B) Position the paddle end (1) of the template on the baby's forehead in the desired position opposite to the cleft side nostril, as shown. (consult referring medical professional for confirmation of placement and optimal angle)

Feed the template tail through the hook loop (2) and slide the tail through until the hook is positioned into the nostril and touching the skin. Fold and crease the extra tail material over and mark it to the desired length (3), as shown by the dotted line (4).

Cut any extra material (5) off of the tail of the template.

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
C) Fold the tail over and trim any excess tail material off to achieve the desired length. (Keep and use this template as a guide for cutting the adhesive tethers to size for future applications during the remainder of the Nasal Elevation System treatment period.)

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A) Lay the paper template over the adhesive tether and cut any excess tether tail material off so the adhesive tether length matches the template length.

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

B) Separate the peel back liner at the pre-cut section, and fold the peel back liner back onto itself as shown so the backing end meets with the cut end of the adhesive tether. (Peeled back liner should be folded in half)

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C) Thread the tail of the adhesive tether through the loop on the nasal hook so the half way point (portion where peel back liner has been folded back) is positioned at the loop opening in the hook.

Pinch the adhesive tether and hook together slightly to stick the hook and tether in this position.

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
D) Peel the remainder of the tail's peel back liner away from the adhesive tether.

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

E) Fold the adhesive back onto itself and pinch the adhesive tether firmly to the hook's loop region. The hook is now secured to the tether.

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A) Fix a skin-barrier tape to the baby's forehead: Based on the position on the forehead of the adhesive tether's paddle, cut a piece of skin barrier larger than the paddle and stick the skin barrier to the forehead. The adhesive tether will be attached to the skin barrier and the skin barrier may stay on the skin for numerous days.

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B) Hold the suspension tape by the paddle-end and gently insert hook into the cleft nostril.

Gently lift until the nostril skin (A) blanches slightly (lightens in color), then determine best suspension angle (B). The most common angle of pull is between vertical and slightly diagonal towards the opposite side of the affected nostril.

When optimal suspension angle has been positioned, and the tension (C) is producing slight blanching of the nostril, remove remaining peel back liner from the paddle end, position the paddle-end over the skin barrier, then press down to adhere the paddle end to the skin barrier.

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C) Rub over the adhered paddle-end for 10-20 seconds to ensure a good adhesion to the skin barrier. The Nasal Elevation installation is complete. Always follow safety precautions to prevent choking on a loose or grasped device.

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